3.7.1 Camera Based Beam Near-Field Propagation Analyzer: M²


3.7.1.1 BeamSquared®

- ISO compliant
- Automatically measure your beam quality in under 1 minutes
- Tune your laser for best operation
- Specifically developed for continuous usage
- Unequaled accuracy using patented Ultracal™ Calibration
- Long optical train & automatic attenuation adjustment
- Flexible mounting configurations, install horizontal or vertically
- Pulsed and CW for most beam diameters and powers
- Compact and portable
- Detectors from 266nm to 10.6µm

The BeamSquared® system is a compact and fully automated tool for measuring the propagation characteristics of CW and pulsed laser systems from the UV to NIR to Telecom wavelengths. Users can also measure wavelengths above 1.8 microns, including CO₂ and terahertz in manual mode (a bench set-up; without the automated optical train) with a Pyrocam™ IV or IIIHR. Our longer optical train and patented Ultracal™ Calibration makes BeamSquared the most accurate product on the market and is ISO 11146 compliant. Its operational robustness and reliability ensures continuous use applications in industry, science, research and development.


Automatic M² - at Production Speeds

The BeamSquared optical train uses a fixed position lens with movable mirrors and camera. The mirrors that direct the focused beam into the camera are moved to precise locations, translating the beam through the near field, the waist, and the far field regions. All these measurements and translations, as well as incremental beam attenuation, are automatically controlled by the BeamSquared software. Design improvements in the BeamSquared system have decreased the measurement reporting time by 2-3 times, making it possible to report M2 in under a minute.



Manual M²

Manual mode is available for wavelengths greater than NIR, particularly Terahertz and above, and for beams that are too large or too small for the BeamSquared optical system. Users are required to provide a manual translation/attenuation apparatus.





Features	
Measurements	
	M2x, M2y, Kx, Ky, BPPx, BPPy
	Width at waist Wx, Wy
	Divergence angle Qx, Qy
	Waist location Zx, Zy
	Rayleigh X, Y
	Astigmatism
	Asymmetry ratio
	Statistical results are available on all measurements
Supports both automated and r	manuai runs
New Hardware	Compays Options includes CD000 Varia DuragenTM III LID or IV
	Camera Options include: SP920, Xeva, Pyrocam™ III HR or IV
	RF Lens Reader
	 Lens must be present for operation Lens configuration data stored with lens (Focal length, calibration wavelength, material, etc.)
	Shutter only open when in live mode
Pupporto hardwara Trigger	Table and attenuator calibration at startup (homing before each run)
Supports hardware Trigger Faster run times than M2-200s	
Vew Interface	
New Interface	Selectable theme colors
	Splash screen with progress bar
2D display	opiasi i scient with progress bai
. В спортау	Selectable Color Palette
	Manual Cursor when not running (Cursor at centroid otherwise)
Caustic Display	Wardar Guisor When the training (Guisor at centrola difference)
Daustic Display	Selecting individual frames
	Auto Aperture
	Exclude points from run
Run Info Display	Excitate points normal
idit ii ilo Diopidy	Displays Caution Notice when beams are non-conforming: (too dark, too bright, misaligned, too large or too small)
	Option to ignore misaligned beams
Editable Settings (Wavelength, I	aser to box distance, Laser to lens and focal length in manual mode)
Calculations	to box distance, East to is in and recall origin in maintain motory
	Frame Results (Total, Min, Peak, % in Aperture, Avg Pwr Density, Beam Width, Centroid, Peak, Cross Sectional Area)
	Laser Results (Waist Width, Divergence, Waist Location Rayleigh Length, M2, K, BPP, Astigmatism, Asymmetry)
	After Lens Results (Waist Width, Divergence, Waist Location Rayleigh Length, Astigmatism, Asymmetry)
	Effective Focal Length of lens
	Fitted/Measured Divergence
	Supported Beam Width calculations
	D4 Sigma
	Knife Edge 10/90 and Programmable
	EPSA - Encircled Power Smallest Aperture (power in a bucket)
Multiple Runs	· · · · · · · · · · · · · · · · · · ·
	Result statistics
	Progress Indicator
Single Page Report	
	Setup information
	Results
	Statistics
	Caustic chart
_ogging/Export data	
	.CVS File



Accuracy by Design

Spiricon products are known for accuracy. Using our patented Ultracal calibration method, auto aperturing to exclude noise beyond the wings of the laser beam, and long optical path, assures the user of the most accurate measurements in the industry.


Designed by Our Customers

Guided by customer input from our widely deployed previous generation M2-200s system, Spiricon redesigned the BeamSquared® to meet the challenging demands of the laser industry. The new BeamSquared system has significantly higher durability and operational robustness for continuous use in a three shifts a day, seven days a week environment. The rigid baseplate and internal optics greatly simplifies and reduces the time for initial set-up and alignment. The lens configuration data is now stored using an RF ID chip embedded in the lens holder which is uploaded automatically by the BeamSquared system when the lens cartridge is inserted in the system, eliminating the need for our customers to keep track of configuration file. Both novice and seasoned users will appreciate these new features along with the time-tested excellence that Spiricon has provided over the years.

Measurements

BeamSquared measures propagation characteristics in both the X and Y axes and displays the following parameters:


- Waist diameters
- Full angle Divergences
- Waist locations
- Rayleigh lengths
- M² or K and BPP factors
- Astigmatism
- Asymmetry



To optimize bench space. BeamSquared can be mounted either horizontally or vertically. Laser beam input port is the same dimension with either mounting method, X = Y, and the same as the M²-200s that it is replacing.

Main Screen Functions

This window displays quantitative measurements of the laser parameters. These include the X and Y beam widths, M2 or K, the divergence angles, the Rayleigh range, and other parameters shown.



This window presents measurements of beam width vs. position for a given run. After measuring a few points, the software extrapolates a curve fit. The Xs and Ys represent individual measurement points. The solid lines present the best fit hyperbola of the beam propagation equation to the measured points. The M² and other laser parameters are computed from the best fit hyperbola since it provides a smoothing of the data points.

This window displays the 2D or 3D beam profile of the currently measured point in the beam propagation curve. This image enables visual intuitive verification of the beam profile behavior through focus. After each run the user can click any individual measured point and observe the beam profile. Outlying or anomalous points can be automatically or manually excluded from the curve fit calculations for more accurate results.



Specifications

Model	BSQ-SP920	BSQ-XC130	BSQ-A	BSQ-PY-M
General				
Camera model	SP920	XC130	BeamSquared® software, software	Pyrocam™ IIIHR or Pyrocam™ IV
Sensor type	Silicon CCD	InGaAs CCD	license, and optical train,	Software only, camera and optical train not included. See individual camera data sheets
Wavelengths	266 – 1100nm	900 – 1700nm	no camera included	
Active area	7.1mm x 5.3mm	9.6mm x 7.6mm		
Elements	1624 x 1224	320 x 256		
Effective pixel	4.4µm x 4.4µm	30μm x 30μm		
Dynamic range	60dB	68dB		
Frame rate	15 fps	100 fps		
Interface	USB 2.0 and 3.0	100 lpc		
Accuracy	±5% typical, ±10% waist location	and Rayleigh length typical		
Measurement cycle time	<1 minute typical, depending on se			
Camera attachment	Standard C-mount, 90° camera or			
Translation system	Standard C-modnit, 90 Camera or Step-motor driven ball screw	T AXIS TOTALIOTT		
Resolution	0.05mm			
	CE, China RoHS			
Compliance	CE, China Rons			
Standard optics	000 440 111/500 51	1000 1700	000 440 111/500 51	N1/A
Lenses included (1)	266-440nm UV 500mm FL 430-700nm VIS 500mm FL 430-700nm VIS 400mm FL 650-1000nm NIR 400mm FL 1000-1700nm Extended NIR 400mm FL	1000-1700nm Extended NIR 400mm FL 650-1000nm NIR 400 FL	266-440nm UV 500mm FL 430-700nm VIS 500mm FL 430-700nm VIS 400mm FL 650-1000nm NIR 400mm FL 1000-1700nm Extended NIR 400mm FL	N/A
Attenuation range				
	Nominally from ND 1.0 to ND 4.8.	Actual values vary with wavelength		N/A
Damage limits (2)	,	,		
	0.15 mW/cm ² CW mode 1.0 µJ/cm ² pulse mode Both of the above for an M ² =1 @ 1064nm	100 mW/cm ²	Depends on type of the camera	See camera data sheets
Optical limits				
Wavelength range	266 to 1100nm	900 to 1700nm	Depends on type of the camera	1.06 to 3000µm
Beam size	BeamSquared Auto Mode 1mm – Varies with wavelength, waist size,			Pyrocam IIIHR 0.8mm – 10mm ma Pyrocam IV 0.8mm – 20mm max Depends on customer mechanics and lens
Minimum beam width	44µm	300µm	N/A	800µm
Software		· ·		
BeamSquared Software	Fast scan method (1 minute) for au	utomatic (ISO) and manual M2 mea	surement	
Environmental				
Storage temperature	-30° C to 65° C			N/A
Storage humidity	95% maximum (non-condensing)			N/A
Operating temperature	10° C to 40° C			N/A
Operating humidity	95% maximum (non-condensing)			N/A
Power requirements (3)				
Input voltage	90 – 264 V AC			N/A
AC Line current	1.6 A			N/A
Line frequency	47Hz to 63Hz			N/A
Physical				
Weight	26 lbs. w/o camera			N/A
Dimensions	217.2mm X 459.5mm X 156.3mm			N/A
Ordering information				
Part Number	SP90502	SP90444	SP90445	SP90410
Notes:	Different lenses are required for different (2) CCD cameras can be damaged by pow While it may be that the laser input pow Use caution and error on the side of saf	wavelength regions, spot sizes and diverge er in excess of 0.15 mW/cm² or energy in e er or energy measures well below this dame by. CCD cameras can be costly to repair or later supplies the power for the system comp	ences, Additional lenses must be ordered separ excess of 1 μ J/cm². BeamSquared employs a frage threshold, it can easily exceed these levels r replace.	ately. ocusing optic. when focused onto the camera sensor.



Accessories Ordering Information

Item	Description	P/N
BSQ-SP920-A	An SP920 camera licensed for BeamSquared®. Sold as an accessory for those also purchasing a BSQ-XC130	SP90521
BSQ-Lens Kit 266-1550	Lens kit that includes 5 BeamSquared lenses: 500mm UV, 500mm VIS, 400mm VIS, 400mm NIR, 400mm XNIR	SP90449
BSQ-Lens Kit 650-1700	Lens kit that includes 2 BeamSquared lenses: 400mm NIR, and 400mm XNIR.	SP90450
BSQ-Lens Kit UV 500mm	Single BeamSquared lens, 500mm focal length, A/R coated for 266-440nm	SP90451
BSQ-Lens Kit VIS 500mm	Single BeamSquared lens, 500mm focal length, A/R coated for 430-700nm	SP90452
BSQ-Lens Kit VIS 400mm	Single BeamSquared lens, 400mm focal length, A/R coated for 430-700nm	SP90453
BSQ-Lens Kit NIR 400mm	Single BeamSquared lens, 400mm focal length, A/R coated for 650-1000nm	SP90454
BSQ-Lens Kit Extended NIR 400mm	Single BeamSquared lens, 400mm focal length, A/R coated for 1000-1550nm	SP90455
BSQ-Lens Kit Extended NIR 600mm	Single BeamSquared lens, 600mm focal length, A/R coated for 1000-1550nm	SP90485
BSQ/BGS-KEY	Includes BeamGage Standard software license in addition to BeamSquared software license	SP90507
BSQ/BGP-KEY	Includes BeamGage Professional software license in addition to BeamSquared software license	SP90508
BSQ SP300 to SP920 upgrade	Camera upgrade	SP90511

