

3.6.1 Camera Based Beam Propagation Analyzer: M^2



M^2 -200s

- Automatically measure your beam quality in under 2 minutes
- Tune your laser for best operation
- ISO compliant
- Specifically developed for continuous usage
- Unequaled accuracy using patented Ultracal™
- Calibration
- Automatic attenuation adjustment
- Pulsed and CW for most beam diameters and powers
- Compact and portable

Not all commercial M^2 measuring instruments conform to the ISO 11146 method of employing a fixed position lens and moving detector. Instead, some manufacturers use a fixed position detector and a moving lens. If the laser beam is diverging or converging within the travel range of a moving lens, the reported M^2 value and other results can be significantly compromised. Spiricon's M^2 -200s and M^2 -200 Beam Propagation Analyzers are fully ISO 11146 compliant.


Automatic M^2 - at Production Speeds

The M^2 -200s optical train uses a fixed position lens and camera. The mirrors that direct the focused beam into the camera are moved to precise locations, translating the beam through both the waist region and the far field regions. All these measurements and translations, as well as incremental beam attenuation, are automatically controlled by the M^2 -200s software. Software improvements in the M^2 -200s, including more efficient algorithm execution, has decreased the measurement reporting time by 2-3 times, making it possible to report M^2 in under two minutes.



Manual M²

Manual mode is available for beams that are too large or too small or at wavelengths outside the standard optical train.



Accuracy by Design


Spiricon products are known for accuracy. Using our patented Ultracal™ calibration method and auto aperturing to exclude noise beyond the wings of the laser beam, assures the user of the most accurate measurements in the industry.

Designed by Our Customers

Spiricon has redesigned the M²-200, the world's top selling beam propagation system to include customer input, increased attention to durability, and operational robustness for continuous use applications - three shifts a day, seven days a week. Novice and seasoned users will appreciate these new features along with the time-tested excellence that the Spiricon M²-200 measurement system has provided over the years.

Main Screen Functions

This window displays quantitative measurements of the laser parameters. These include the X and Y beam widths, M² or K, the divergence angles, the Rayleigh range, and other parameters shown.



This window presents measurements of beam width vs. position for a given run. After measuring a few points, the software extrapolates a curve fit. The Xs and Ys represent individual measurement points. The solid lines present the best fit hyperbola of the beam propagation equation to the measured points. The M² and other laser parameters are computed from the best fit hyperbola since it provides a smoothing of the data points.

The 2D or 3D beam profile of the currently measured point in the beam propagation curve. This image enables visual intuitive verification of the beam profile behavior through focus. After each run the user can click any individual measured point and observe the beam profile. Outlying or anomalous points can be automatically or manually excluded from the curve fit calculations for more accurate results.

3.6.1.1 Specifications for the M²-200s

General		
Accuracy	±5% typical, ±12% waist location and Rayleigh length typical (Note: Accuracy can be degraded by a variety of situations)	
Measurement Cycle Time	2-3 minutes typical, depending on setup conditions and operating mode	
Camera Attachment	Std C-mount, 90° camera on axis rotation	
Translation System	Step motor-driven lead screw	
Translation Pitch	4 mm/rev optical pitch	
Step Angle	1.8° (200 steps/rev)	
Sample Range	M ² - 200 s 190 - 600 mm, typical	
Camera Specifications (for GRAS20 camera)		
Imager	1/1.8" CCD, 1600 x 1200 pixels	
Dynamic Range	12 bit A to D	
Frame Rates	7.5 FPS (at full resolution)	
Pixel size	4.4µm x 4.4µm	
Gain	0 to 25dB	
Shutter Control	Programmable from 110µs to 70ms	
S/N Ratio	59dB at min gain	
Trigger Input	Edge sensitive 3.3 / 5Vdc LVTTTL / TTL (positive or negative, user programmable)	
	Minimum pulse width 10µs	
Trigger Out	External Trigger cable provided	
Voltage Requirement	3.3Vdc LVTTTL, Programmable	
Power Consumption	Powered over Firewire Cable <3.5watts	
Dimensions	44mm (1.74") wide, 29mm (1.14") tall and 66mm(2.6") deep	
Mass	104g (3.7oz)	
Environmental		
Storage Temperature	-30°C to 65°C	
Storage Humidity	95% maximum (non-condensing)	
Operating Temperature	10°C to 40°C	
Operating Humidity	95% maximum (non-condensing)	
Power Requirements*		
Line Voltage	95V AC to 250V AC	
Line Frequency	47Hz to 63Hz	
Maximum Power	4.5 Watts	
* For the Optical Train only. The PC computer supplies the power for the system components, such as the CCD camera. An external power supply is provided for Laptop computer use.		
Physical		
Weight	M ² -200s... 6.8 kg (without camera)	
Measurements		
(Statistical results are available on all measurements)	M ² x, M ² y, Kx, Ky, BPPx, BPPy Width at waist Wx, Wy Divergence angle qx, qy Waist location Zx, Zy Rayleigh X, Y Astigmatism Asymmetry ratio	
Wavelength Range		
Different lenses are needed for different wavelength regions The M ² -200s model include 3 standard lenses with nominal 300mm focal lengths. See below		
M ² -200s-FW	266 - 587nm (included) 400 - 750nm (included) 650 - 1300nm (included) 1000 - 1300nm (optional)	
Attenuation Range Nominally from ND 0 to ND 4.8. Actual values vary with wavelength		
Beam Size	0.5mm - 10mm M ² -200s Varies with wavelength, waist size and location, and M ²	
Damage Limits ¹		
Camera	0.15 µW/cm ² CW mode for a 10 mm input beam diameter 1.0 µJ/cm ² pulse mode for a 10 mm input beam diameter Both of the above for an M ² = 1 @ 1064nm	
¹ CCD cameras can be damaged by power in excess of 100 mW/cm ² or energy in excess of 100 mJ/cm ² . The M ² -200s employs a focusing optic. While it may be that the laser input power or energy measures well below this damage threshold, it can easily exceed these levels when focused onto the camera sensor. Use caution and error on the side of safety. CCD cameras can be costly to repair or replace.		
Ordering Information		
Item	Description	P/N
M²-200s Beam Propagation Analyzer		
M ² -200s-FW	M ² -200 software, software license, GRAS 20 Firewire camera, short optical train, automatic and manual operation, recommended for 266nm - 1064nm wavelengths	SP90144
M ² -200s-FW-A	M ² -200 software, software license, short optical train, automatic and manual operation, recommended for 266nm - 1064nm wavelengths (GRAS 20 camera not included)	SP90145
M ² -200sM-FW	Manual mode M ² -200 software, software license, GRAS 20 Firewire camera, manual operation with a GRAS 20 Firewire camera (optical train not included)	SP90146
M ² -200sM-FW-A	Manual mode M ² -200 software, software license, manual operation with a Firewire camera (GRAS 20 Firewire camera and optical train not included)	SP90147